


# EXINORE Sleepfit



## SPECIFICALLY INTENDED FOR

- Anyone who is irritated by snoring
- People who work night shifts
- Students who want to be able to concentrate better
- Anyone who finds themselves in a noisy environment
- Travelling by plane, car, bus or train

EXINORE SleepFit has been given a face-lift. With its new shape, this custom-made earpiece now fits even better in the auditory canal giving optimal wearer comfort due to use of flexible material. The SleepFit provides a perfect fit which eliminates unwanted sound from your environment. Although noise is suppressed significantly, the SleepFit still allows a functional sound level (also see: our advice). Therefore an alarm clock or fire alarms can still be heard while wearing the SleepFit.



The SleepFit offers a benefit in many different situations. Students who particularly need to concentrate use them. People who work night shifts and want to be able to sleep during the day without being disturbed. This also applies to people who live close to entertainment centers. The SleepFit is also ideal for those who travel frequently for business and other purposes. This custom-made product provides the rest you need, when you need it.

### Our advice:

There are people who think they will not be able to hear anything if they wear the SleepFit. That is a persistent misconception. Because sounds enter the body via the nose, mouth and bones as well (bone conduction), some sounds can still be perceived. The noise of an alarm clock, for instance, or a fire alarm will still be effective. By clarifying this to your customers, you can remove their concern.

### Our advice:

Some people might need to get used to wearing earplugs. Adapting can take a week, but doesn't live up to the benefits it will provide for the years after that.

## COLOUR

The SleepFit is available in either transparent or yellow.

**Note! When no colour has been indicated the SleepFit will automatically be provided in yellow. No other colours are possible due to the choice of material which guarantees the softness of the SleepFit.**

## MATERIAL

The SleepFit is made from extra soft silicone material, with a hardness degree of 40 Shore.

## BENEFITS

- Custom-made, so it fits perfectly
- No protruding parts; fits entirely inside the auditory canal
- Made from extra soft, light-weight material
- Easy to insert and remove
- Comfortable
- Quick and easy to clean
- Improves the quality of sleep
- Better concentration due to attenuation of surrounding noises
- Alarm signals from an alarm clock or fire alarm can still be heard



**EXINORE**  
Hearing protection

[www.exinore.com](http://www.exinore.com)